

Transmission electron microscopy (TEM) enables characterization of powders and thin films (which can be prepared in a target preparation from bulk materials) by direct imaging with up to atomic resolution. The image information can be locally correlated with spectroscopic techniques (EELS/EFTEM and EDX) to provide semi-quantitative elemental composition/maps with sub-nanometer resolution. All of these techniques can also be performed in-situ, e.g. during heating, electrical biasing or straining to directly correlate structural changes and materials properties. For complex three-dimensional structures, electron tomography can be used to generate a 3D representation of the material with a spatial resolution of 1–2 nm.

Contact

Dr. Christian Kübel

Phone +49 721 608-28970, fax +49 721 608-28976, email christian.kuebel@kit.edu

Institute of Nanotechnology (INT) - www.int.kit.edu/english

Features

- FEI Titan 80–300 (aberration corrected TEM)
- Resolution:
 - 0.08 nm information limit TEM
 - 0.14 nm resolution in STEM
 - 0.7 eV energy resolution EELS
- Imaging and Analysis Techniques:
 - BF-TEM, aberration corrected HRTEM
 - HAADF-STEM, HRSTEM
 - EFTEM, EELS, EDX
 - (S)TEM tomography
 - electron diffraction, electron precession
 - orientation mapping
 - Lorentz imaging
 - low-dose techniques & cryo imaging
- In-situ Techniques:
 - Heating (Protochips Aduro: RT-1200 °C; Gatan 652: RT-800 °C)
 - Cooling (Gatan 915: LN₂-80 °C)
 - Straining (Hysitron Picoindenter PI 95 and Gatan 654)
 - Electrical Biasing (Protochips Aduro)
 - Electro chemistry (Protochips Poseidon 500)
- Sample preparation:
 - Thin films or nano powders can be directly analyzed without additional preparation
 - Target preparation by FIB lift-out (for details see FIB description) with final polishing by low-voltage Argon ion beam (Fischione 1040 NanoMill)
 - Electro polishing
 - Classical preparation by cutting, grinding, argon ion milling or microtomy

Limitations/constraints

- Sample has to be a solid at LN₂ temperatures and stable under high vacuum conditions
- Maximum sample thickness: 10–2000 nm (depending on resolution and technique)
- Depending on the structure and chemical composition, the sample might be sensitive to the electron beam resulting in changes during analysis
- Except in tomography, TEM always provides an image/analysis of the projected structure of a sample
- H, He und Li cannot be detected by our analytical techniques

Typical results




Fig. 1: HAADF-STEM image (filtered by NAD) of a $La_{1-x}Sr_xMnO_3/SrTiO_3$ interface with the individual atomic columns well resolved across the interface. Overlaid is an EELS/EDX intensity profile across this interface. P.M. Leufke and D. Wang et al., *Thin solid films*, 2012, 520, 5521-5527.




Fig. 3: Geometric phase analysis reveals the local strain distribution around the triple line in the image above. H. Rösner and C. Kübel et al., *Acta Mat.*, 2011, 59, 7380-7387.




Fig. 2: Atomic resolution TEM image of a triple and a quadruple line at the interface between $\Sigma 3$ boundaries and a $\Sigma 9$ boundary in nanocrystalline palladium. H. Rösner and C. Kübel et al., *Acta Mat.*, 2011, 59, 7380-7387.




Fig. 4: In-situ orientation mapping (different color correspond to different crystal orientations) of the grain structure changes in nanocrystalline gold during straining – selected images of the straining series showing anomalous grain growth. A. Kobler and C. Kübel et al., *Ultramicroscopy*, 2013, 128, 68-81.

Typical results (*continued*)




Fig. 5: EFTEM mapping (Si-blue, C-red) and HRTEM image of nanocrystalline silicon particles with a covalently bound C18 shell. The EFTEM maps reveal the ~ 1.2 nm wide carbon shell around the silicon core. Sample provided by G. Ozin, University of Toronto.




Fig. 7: HRTEM image of nano graphene with the corresponding low-loss EELS spectrum showing the characteristic π and $\pi+\sigma$ plasmon losses. J. Biener and D. Wang et al., *Adv. Mater.* 2012, 24, 5083–5087.




Fig. 6: HAADF-STEM image with EDX compositional mapping of the different layers in a silicon quantum dot based organic LED (SiLED). F. Maier-Flaig and C. Kübel et al., *Nano Letters*, 2013, online; DOI: 10.1021/nl400975u.




Fig. 8: HRTEM image of a Fe/LiF/C anode for lithium ion batteries revealing α -iron nanoparticles each surrounded by a few graphene layers. R. Prakash and C. Kübel et al., *J. Power Sources*, 2011, 196, 5936–5944.

Typical results (*continued*)




Fig. 9: Electron tomographic reconstruction of a self-assembled CdS nano cluster superlattice (with additional 5 nm gold particles in yellow). The two digital slices, one unit cell apart, show a single vacancy, an extended vacancy and dislocations in 3D. T. Levchenko and C. Kübel et al., *Chem. Eur. J.*, 2011, 17, 14394-14398.




Fig. 10: HAADF-STEM and HRTEM imaging of uniform ThO₂ nanorods. D. Hudry and E. Courtois et al., *Chem. Eur. J.*, 2013, 19(17), 5297–5305.